

The 4000 curved stairlift

Safe, reliable and easy to use. Custom-made to fit perfectly in your home for a neat finish, thanks to our advanced digital staircase measurement system. Easy to use and keep clean, as well as to maintain and service. Smooth ride, strength and stability thanks to state-of-the-art twin-rail design technology.

Perfect for compact stairs and hallways, and even spiral staircases.



Handicare is proud to have achieved ISO certifications for quality, environment and health and safety. We are committed to providing stairlifts which are easy to use, reliable, comfortable and which meet the highest possible quality and safety standards. All new Handicare stairlifts are delivered with a comprehensive manufacturer's warranty and are independently tested and approved by Europe's leading certification organisation including Liftinstituut.

www.handicare-stairlifts.co.uk

Style seat White Beige Red	Smart seat Sand Cocoa Ruby Sate	<image/>	Simplicity seat*
Rail Colours	Light Grey Jet Black Brown	TAILORED: Eas powered options today and tomo	sily upgrade with s to suit your needs rrow
SPECIFICATIONS		COMPACT: Folds to leave more room ↓ ∠ on the stairs when not in use. Hinged ∧ K rail allows clear access on the lower	
Applications	Residential indoor and outdoor		
Capacity	Standard 136kg (21½st) Heavy Duty 160kg (25½st)	SAFETY: Retractable seat belt, emergency stop, and safety sensors	
		STYLISH: Modern sleek design with a choice of seats and colours	
Minimum track intrusion into staircase	125mm	OUTDOOR: With improved UV protection and a water-resistant cover (IP44 rated)	
Speed	0.12m/sec		
Maximum staircase angle	56°	B C costs and batter c if there's a powe	CIENT: Low running ries to remain operational er cut
*Simplicity seat is supplied as standard			

*Simplicity seat is supplied as standard

The certified Handicare partner near you:



Please be aware that actual colours might vary slightly from those shown. All information is provided to you in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.