


M-Series Ceiling Lift Model Comparison

| |  |  |  |
|--|--|---|---|
| | Essential | Clinical | Bariatric |
| Lithium-ion battery technology | ✓ | ✓ | ✓ |
| SWL options: 130, 200 or 272 kg (286, 440 or 600 lbs) | ✓ | ✓ | |
| SWL options: 363, 454 or 544 kg (800, 1000 or 1200 lbs) | | | ✓ |
| Manual traverse option | ✓ | ✓ | |
| Powered traverse option | ✓ | ✓ | ✓ |
| In-track charging option | ✓ | ✓ | ✓ |
| Models available for competitors' track, including Arjo, Guldmann and Liko | ✓ | ✓ | ✓ |
| 3-year warranty including battery | ✓ | ✓ | ✓ |
| Return to Charge (RTC) & Return to Home (RTH) options | ✓ | ✓ | |
| Digital display and menu: choice of on lift, or on handset | | ✓ | ✓ |
| Cycle counter and Advanced Life Analytics | | ✓ | ✓ |
| "3-Down" three-mode safety function | | ✓ | ✓ |
| Bluetooth/Wi-Fi Connectivity option | | ✓ | ✓ |
| 360-degree rotation | | | ✓ |
| Motorized 4-point tilting carry bar option | | | ✓ |